




TRIOPTICS

TRIOPTICS GMBH · OPTISCHE INSTRUMENTE

OptiCentric®

Alignment and Cementing
Workstation



APPLICATIONS

- Alignment and cementing of lenses on an arbor as a preliminary step for lens edging
- Centration error measurement of single lenses and achromats (optional)

FEATURES

- Centration error measurement in submicron range
- Graphical and numerical display of the alignment path
- Option: Fully automated alignment of the lens to the reference axis with the help of the SmartAlign® Algorithm and piezo manipulators
- Flexible system adaptable to different production environments



TRIOPTICS